
Rheomuco

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
- ・ Serisier, D.J., Carroll, M.P., Shute, J.K. et al. Macrorheology of cystic fibrosis, chronic obstructive pulmonary disease & normal sputum. *Respir Res* 10, 63 (2009).
- ・ Hurst, J.R., Wedzicha, J.A. Management and prevention of chronic obstructive pulmonary disease exacerbations: a state of the art review. *BMC Med* 7, 40 (2009).
- ・ Nettle, C., Jenkins, L., Curtis, D., Badiei, N., Lewis, K., Williams, P. & Daniels, D. (2017). Linear rheology as a potential monitoring tool for sputum in patients with Chronic Obstructive Pulmonary Disease (COPD). *Biorheology*, 1-14.

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Rheomuco by Rheonova CONFIDENTIAL DID\MEARP\V1.0

Rheomuco Measurement report

Measurement parameters		Device specifications	
Date - time	10 Jun 2020 - 4:16 PM	Device serial number	501
Measurement identification	test EN	Motors firmware	1.34
Operator (profile)	PB (Administrator)	Coders firmware	1.17
Temperature setpoint	37 °C	Last maintenance date	19 May 2019
Gap	1 mm	Operating software version	1.4.0.0 (medical)
Geometries	rough, 25 mm	Measuring procedure version	Strain sweep n.1
Reinstalled sample	no		

Measurement result		Value
Measurement validity ⁽¹⁾	Valid results: linear region picked at 0.05 strain, flow point found	
Viscoelastic modulus	G^*	3.32 Pa
Damping factor	$\tan(\delta)$	0.338
Critical strain	γ_c	7.5
Critical stress	σ_c	16.6 Pa

⁽¹⁾ The measurement is valid if and only if no software or hardware incoherence are detected by the operating software during the analysis or by the algorithm during results treatment (see page 2)

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